

Fall Management

Falls are a well-recognized problem within nursing homes and assisted living settings. Each year, thousands of residents fall, which can lead to serious injuries and even death. Monitor resident movement and prevent falls with PalatiumCare and integrated Fall Management devices.

PalatiumCare offers three pressure sensitive pads that integrate with your wire-less nurse call system: Chair Pad, Bed Pad, and Floor Pad. The Chair and Bed Pads alert caregivers when a resident is moving off the pad, and the Floor Pad alerts as soon as a resident steps on to the mat.

Easily integrate each pad with your PalatiumCare nurse call system:

For more information on PalatiumCare and integrated Fall Management devices, call 1-888-725-2848. See more PalatiumCare resident safety and call devices at www.PalatiumCare.com


Select the PalatiumCare integrated fall management device that best meets your residents' needs.

Call 888-725-2848 to consult a sales representative.

Chair Pad



Bed Pad



Floor Pad



Connect the pad to the fall management monitor using the Sensor Pad port.

Connect the universal transmitter to the fall management monitor using the Nurse Call port.

Use the Reset button on the fall management monitor to learn the pad, monitor and transmitter into your PalatiumCare system.

Disable the audible alert to eliminate nuisance alerts.