

# PALATIUM CARE

## Crest Fall Fighter® | Bed Pad

Falls are a well recognized problem within nursing homes and assisted living settings. Each year, thousands of residents experience falls, which can lead to serious injuries or even death. After a fall occurs, senior living communities face additional costs associated with the extra time and care required with a resident who has a serious injury. Avoid these extra costs and stress on your staff and residents by monitoring resident movement with a Crest Fall Fighter® Bed Pad and PalatiumCare Universal Transmitter (UTX).



Connect the Bed Pad to the Crest Fall Fighter® Monitor using the Sensor Pad port

Connect the UTX to the Crest Fall Fighter® Monitor using the Nurse Call port



When the UTX is learned into the PalatiumCare system, the Bed Pad will alert staff through your nurse call system. Disable the audible alert on the Fall Fighter® Monitor to eliminate nuisance alarms.

For more information on PalatiumCare and integrated Fall Fighter® devices, call 1-888-725-2848. Additional PalatiumCare resident safety and call devices can be viewed at [www.PalatiumCare.com](http://www.PalatiumCare.com).

